



Instructions for Use

Did you know you can keep your dressings or PICC lines dry while showering or bathing? You can with Seal-Tight Protector... *The easy durable way to stay dry.*

Here's how easy it is!

1. After removing protector from the box, unfold completely.
2. Pull SEAL-TIGHT® Protector over dressing or PICC line by stretching the opening of the protector over your dressing while pulling up. You will know it's at the right position on your arm or leg when the seal opening sits comfortably above the top edge of the dressing.

Important note: The seal should be tight, after all, it's keeping water off of your dressing. Never trim with scissors or puncture. It's ok if the bag appears oversized, as this allows room for large dressings and easy movement.

3. Partially submerge protector into the water while simultaneously holding a portion of the proximal (uppermost) seal away from your skin with one finger. This will allow any excess air to escape easily from the bag.
4. When excess air has escaped, remove your finger so the seal fits snugly against your skin.
5. To Remove: Dry the protector thoroughly before stretching the seal over dressing and remove. **DO NOT** remove by pulling on the bag.

For therapeutic care/daily showering or bathing. Use care around water, may be slippery when wet. Consult your doctor to be sure swimming is appropriate in your situation.

Brownmed warrants this product to be free of manufacturing defects. Intended for use on single patient.

