

Multi-Pads
Provide an anatomical support pad that can be used as medial/lateral posting orthotic or neutral heel cup. Multi-Pad can be modified through use of accessory Posting Wedges to provide varus and valgus buildup. Available in sizes Small, Medium and Large, with 2 per package.



Posting Wedges
May be used with other Soft Stride™ products. Included with Multi-Pads. Maintain original placement due to self-adhering properties. Layer to achieve desired posting. Available in sizes Small, Medium and Large, with 6 per package.



Tender Foot Pads
Separate metatarsal shafts and cushion the metatarsal heads. Unique design alleviates pain and discomfort in forefoot. Aid in treatment of neuromas and metatarsalgia. Available in sizes Small, Medium and Large, with 2 per package.

Thin Insoles
Reduce shock over entire plantar surface with lightweight and thinner forefoot design. Contain extra polymer in heel and arch areas to reduce shock. Available in sizes A, B, C, D or E, with 2 per package.



Underliners
Designed with equal thickness over entire surface, absorbing shock to plantar surface of foot. Relieve pain due to degenerative fat pads. Excellent for walking or working on hard surfaces. Available in sizes A, B, C, D, E, F, G, H, I or J, with 2 per package.

Sheet Products
Prefabricated of viscoelastic polymer, Soft Stride™ Sheets provide an alternative to preformed foot orthotics. Useful for customization or when making an orthotic is desired. Diamond design allows for polymer flow when pressure is applied creating a softer feel. Use for prosthetic applications to provide contact cushioning against stumps, for creating customized orthotics and for shock absorbing mats where needed. Size and thickness variations available, with or without top cover.

Sheet Products: Inquire about size and thickness variations, with or without top covers.

Corporate:
1300 Lundberg Drive West
Spirit Lake, IA 51360 USA
TF 800-843-4395
F 712-336-2874
P +1-712-336-4395

Medical Professional & Retail:
4435 Main Street, Suite 820
Kansas City, MO 64111
TF 877-853-5518
F 816-581-7005



72099

info@brownmed.com • www.brownmed.com

Soft Stride™ Orthotics

Viscoelastic Polymer Products

Soft Stride™ is a comprehensive line of flexible foot orthotics that addresses many specific ailments and injuries to feet, heels, ankles, arches and legs. These products help patients avoid or delay surgery, attenuate shock and provide pain relief. Athletes receive comfort during strenuous training and people who stand for long periods benefit with decreased fatigue and pain.

The viscoelastic polymer mimics the qualities of soft body tissue. This unique polymer provides redistribution of pressure and attenuates impact.

100% Viscoelastic Polymer
- contains no fillers

Impact Energy Control
- reduces fatigue and treats injuries to the musculoskeletal system.

Modular System
- allows healthcare professional to build custom orthotics using biomechanically designed orthotic components.

Non-cellular Material
- prevents bottoming out, compression set and shock wave to skeletal system.

Flexible Cloth Covering
- facilitates impact dispersion and aids in evaporative cooling.

Natural Surface Tack
- features self-adhering bottom that maintains original placement.



Bone Spur Pads

Contain existing soft tissue around calcaneus and relieve direct pressure in heel spur, bruised and sensitive area. Available in sizes Small, Medium and Large, with 2 per package.

Extended Heel Cushions
Absorb heel strike shock and provide complete mid-foot to posterior foot cushion. Are useful in treating plantar fasciitis, heel strike impact reduction and Achilles tendonitis. Available in sizes Small, Medium and Large, with 2 per package.



Heel Wedges

Provide 15 degree declination angle. Useful for leg length discrepancy, Achilles tendonitis and heel strike impact reduction. Available in sizes Small, Medium and Large, with 2 per package.

Longitudinal Arch Supports
Support and maintain normal longitudinal arch. May be used alone in shoe or with Soft Stride™ Thin Insoles. Can be ordered with or without top covering. Order without top cover if used with Thin Insoles. Available in sizes Small, Medium and Large, with 2 per package.



Metatarsal Pads

Support transverse arch and cushion metatarsal heads. May be used with Soft Stride™ Underliners. Can be ordered with or without top covering. Order without top cover if used with Underliners. Available in sizes Small, Medium and Large, with 2 per package.

