

Polar Ice®

Cold Therapy Support

LONG LASTING • EASY TO USE

Encapsulated Ice Pack

- Conforms well to the body
- Transports easily in cooler with ice
- Provides quick refreeze time

Soft Fleece Material

- Eliminates migration of ice wrap during use
- Provides comfortable compression
- Insulates skin to prevent frost bite
- Wicks moisture away



10303 (Small)
10304 (Medium)
10305 (Large)
Finger Sleeve



30011 (Small) **30013** (Large)
30012 (Med) **30014** (X-Large)
Shin Wrap **30036** (Universal)



30100
Shoulder/Hip Wrap



30101
Wrist/Elbow Wrap



30102
CPM Knee Wrap



30103 (Standard)
30104 (Large)
Arthroscopy Knee Wrap



30105
Foot/Ankle Wrap



30107
Quick Wrap



30108
Back Wrap



30109
TMJ Wrap

Use For

- Post-Op Cryotherapy
- Sports Injuries
- Rehabilitation
- Soft Tissue Trauma
- Edema and Swelling

Convenient To Use In

- Hospital
- Home
- Clinic
- Industry
- Sports

ORDER #	ITEM	SIZE
10303	Finger Sleeve	Small
10304	Finger Sleeve	Medium
10305	Finger Sleeve	Large
30011	Shin Wrap - Small	12-14"
30012	Shin Wrap - Medium	14-15"
30013	Shin Wrap - Large	15-17"
30014	Shin Wrap - XLarge	17-19"
30036	Shin Wrap	Universal
30100	Shoulder/Hip Wrap	Universal
30101	Wrist/Elbow Wrap	Universal
30102	CPM Knee Wrap	Universal
30103	Arthroscopy Knee Wrap - Standard	13" x 15.75"
30104	Arthroscopy Knee Wrap - Large	13" x 26.75"
30105	Foot/Ankle Wrap	Universal
30107	Quick Wrap	9" x 12"
30108	Back Wrap	Universal
30109	TMJ Wrap	Universal



Executive:
280 Summer Street
Suite 400
Boston, MA 02210
P 857-317-3354
F 857-317-3355

Sales:
4435 Main Street
Suite 820
Kansas City, MO 64111
TF 877-853-5518
P 816-581-7001
F 816-581-7005

Manufacturing:
1300 Lundberg Drive West
Spirit Lake, IA 51360
TF 800-843-4395
P 712-336-4395
F 712-336-2874

www.brownmed.com

