

Toe Hold®



Adjustable Stabilization of First Metatarsal-Phalangeal Joint

- maintains the hallux in any corrective plane the physician desires (neutral, adducted, varus, dorsiflexed, or plantarflexed position) with appropriate placement of hook and loop strap
- improves hallux purchase during ambulation by stretching soft tissue contractures
- provides support for turf toe or broken hallux
- non-latex foam construction
- one piece design applies easily and stretches to fit without compromising stability

Patent No. D341,424

Effective Treatment Following

- first MPJ realignment
- hallus valgus procedure

ORDER #	SIZE
50201	S/M Right
50202	S/M Left
50203	M/L Right
50204	M/L Left

S/M for shoe widths AA – B
M/L for shoe widths C – EEEE



www.brownmed.com



Executive:

280 Summer Street, Suite 400
Boston, MA 02210
P 857-317-3354
F 857-317-3355

Sales:

4435 Main Street, Suite 820
Kansas City, MO 64111
TF 877-853-5518
P 816-581-7001
F 816-581-7005

Manufacturing:

1300 Lundberg Drive West
Spirit Lake, IA 51360
TF 800-843-4395
P 712-336-4395
F 712-336-2874

50299 Rev 2