

Steady Step[®] Walker



Steady Step[®] Walker is recommended for immobilization and wound care treatment.

- lightweight glass-filled nylon shell
- provides superior protection to injured limb with shell design
- low profile design
- allows convenient and consistent step-in application for medical examinations and patient use
- adaptable for diabetic wound care
- suggested reimbursement code L4386

Patent No. 6,361,514

Standard

ORDER #	SIZE	MEN	WOMEN
52021	Small	3-5½	4-6½
52022	Medium	6-8½	7-9½
52023	Large	9-12	10-13

Low-Top

ORDER #	SIZE	MEN	WOMEN
52031	Small	3-5½	4-6½
52032	Medium	6-8½	7-9½
52033	Large	9-12	10-13

Easily customize initial fitting with included orthotics

- arch support
- removable insole
- shock absorbing pad
- heel stabilizer



www.brownmed.com



Executive:

280 Summer Street
Suite 400
Boston, MA 02210
P 857-317-3354
F 857-317-3355

Sales:

4435 Main Street
Suite 820
Kansas City, MO 64111
TF 877-853-5518
P 816-581-7001
F 816-581-7005

Manufacturing:

1300 Lundberg Drive West
Spirit Lake, IA 51360
P 712-336-4395
F 712-336-2874