



## Viscoelastic Polymer Products

SoftStride™ is a comprehensive line of flexible foot orthotics that addresses many specific ailments and injuries to feet, heels, ankles, arches and legs. These products help patients avoid or delay surgery, attenuate shock and provide pain relief. Athletes receive comfort during strenuous training and people who stand for long periods benefit from decreased fatigue and pain. The viscoelastic polymer mimics the qualities of soft body tissue. This unique polymer provides redistribution of pressure and attenuates impact.



fits any style shoe



unique shock absorbing material



cushions & supports



will not bottom out



natural gripping surface holds pad in place




may be trimmed with scissors to provide better fit

## HINDFOOT PRODUCTS



### Bone Spur Pads

Contain existing soft tissue around calcaneus and relieve direct pressure in heel spur, bruised and sensitive area. Available in sizes Small, Medium and Large, with 2 per package.



### Extended Heel Cushions

Absorb heel strike shock and provide complete mid-foot to posterior foot cushion. Are useful in treating plantar fasciitis, heel strike impact reduction and Achilles tendonitis. Available in sizes Small, Medium and Large, with 2 per package.



### Heel Wedges

Provide 15 degree declination angle. Useful for leg length discrepancy, Achilles tendonitis and heel strike impact reduction. Available in sizes Small, Medium and Large, with 2 per package.



### Longitudinal Arch Supports

Support and maintain normal longitudinal arch. May be used alone in shoe or with SoftStride™ Thin Insoles. Can be ordered with or without top covering. Order without top cover if used with Thin Insoles. Available in sizes Small, Medium and Large, with 2 per package.



### Multi-Pads

Provide an anatomical support pad that can be used as medial/lateral posting orthotic or neutral heel cup. Multi-Pad can be modified through use of accessory Posting Wedges to provide valgus buildup. Available in sizes Small, Medium and Large, with 2 per package.



### Posting Wedges

May be used with other SoftStride™ products. Included with Multi-Pads. Maintain original placement due to self-adhering properties. Layer to achieve desired posting. Available in sizes Small, Medium and Large, with 6 per package.



## SoftStride™ Size Charts

### Bone Spur Pads

Order No.	Size	Men	Women	EUR
71301	Small	4 - 6	5 - 7	36 - 38.5
71302	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71303	Large	11+	12+	44+

### Extended Heel Cushion

Order No.	Size	Men	Women	EUR
71311	Small	4 - 6	5 - 7	36 - 38.5
71312	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71313	Large	11+	12+	44+

### Heel Wedge

Order No.	Size	Men	Women	EUR
71341	Small	4 - 6	5 - 7	36 - 38.5
71342	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71343	Large	11+	12+	44+

### Longitudinal Arch Support

Order No.	Size	Men	Women	EUR
71351	Small	4 - 6	5 - 7	36 - 38.5
71352	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71353	Large	11+	12+	44+

### Multi Pad

Order No.	Size	Men	Women	EUR
71371	Small	4 - 6	5 - 7	36 - 38.5
71372	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71373	Large	11+	12+	44+

### Posting Wedge

Order No.	Size	Men	Women	EUR
71381	Small	4 - 6	5 - 7	36 - 38.5
71382	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71383	Large	11+	12+	44+

### Metatarsal Pads

Order No.	Size	Men	Women	EUR
71361	Small	4 - 6	5 - 7	36 - 38.5
71362	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71363	Large	11+	12+	44+

### Tenderfoot Pad

Order No.	Size	Men	Women	EUR
71411	Small	4 - 6	5 - 7	36 - 38.5
71412	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71413	Large	11+	12+	44+

### Thin Insoles

Order No.	Size	Men	Women	EUR
71421	A	3 - 5	4 - 6	34.5 - 37
71422	B	6 - 8	7 - 9	38.5 - 40.5
71423	C	9 - 11	10 - 12	41.5 - 44
71424	D	12 - 14		45.5 - 48
71425	E	15 - 17		48+


### Underliners

Order No.	Size	Men	Women	EUR
71430	A		4½ - 5½	35 - 36.5
71431	B		6 - 7	37 - 38.5
71432	C	6½ - 7½	7½ - 8½	39 - 40
71433	D	8 - 9	9 - 10	40.5 - 42
71434	E	9½ - 10½	10½ - 11½	42.5 - 43.5
71435	F	11 - 12	12+	44 - 45.5
71436	G	12½ - 13½		46 - 47
71437	H	14 - 15		48 - 49
71438	I	15½ - 16½		50 - 51
71439	J	17 - 18		52+

### Sheet Product

Order No.	Dimension	Thickness
71406	14 x 9.5" / 35.5 x 24.1cm	3/16"

## FOREFOOT PRODUCTS



### Metatarsal Pads

Support transverse arch and cushion metatarsal heads. May be used with SoftStride™ Underliners. Can be ordered with or without top covering. Order without top cover if used with Underliners. Available in sizes Small, Medium and Large, 2 with 2 per package.



### Tender Foot Pads

Separate metatarsal shafts and cushion the metatarsal heads. Unique design alleviates pain and discomfort in forefoot. Aid in treatment of neuromas and metatarsalgia. Available in sizes Small, Medium and Large, with 2 per package.



## INSOLE/UNDERLINDER/SHEET PRODUCTS



### Thin Insoles

Reduce shock over entire plantar surface with lightweight and thinner forefoot design. Contain extra polymer in heel and arch areas to reduce shock. Available in sizes A, B, C, D or E, with 2 per package.



### Underliners

Designed with equal thickness over entire surface, absorbing shock to plantar surface of foot. Relieve pain due to degenerative fat pads. Excellent for walking or working on hard surfaces. Available in sizes A, B, C, D, E, F, G, H, I or J, with 2 per package.



### Sheet Products

Prefabricated of viscoelastic polymer, SoftStride™ Sheets provide an alternative to preformed foot orthotics. Useful for customization or when making an orthotic is desired. Diamond design allows for polymer flow when pressure is applied creating a softer feel. Use for prosthetic applications to provide contact cushioning against stumps, for creating customized orthotics and for shock absorbing mats where needed.



www.brownmed.com



### Executive:

101 Federal Street  
29th Floor  
Boston, MA 02110  
P 857-317-3354  
F 857-317-3355

### Sales:

4435 Main Street, Suite 820  
Kansas City, MO 64111  
TF 877-853-5518  
P 816-581-7001  
F 816-581-7005

### Manufacturing:

1300 Lundberg Drive West  
Spirit Lake, IA 51360  
P 712-336-4395  
F 712-336-2874