



Backed by our 1 million step guarantee.

Viscoelastic Polymer Products

SoftStride™ is a comprehensive line of flexible foot orthotics that addresses many specific ailments and injuries to feet, heels, ankles, arches and legs. These products help patients avoid or delay surgery, attenuate shock and provide pain relief. Athletes receive comfort during strenuous training and people who stand for long periods benefit from decreased fatigue and pain. The viscoelastic polymer mimics the qualities of soft body tissue. This unique polymer provides redistribution of pressure and attenuates impact.



fits in most shoes



unique shock absorbing material



cushions & supports



will not bottom out



position in bottom of shoes, natural gripping surface prevents slipping



may be trimmed with scissors to provide better fit

HINDFOOT PRODUCTS



Bone Spur Pads - Cobra style

Contain existing soft tissue around calcaneus and relieve direct pressure in heel spur, bruised and sensitive area. Available in sizes Small, Medium and Large, with 2 per package.




Bone Spur Pads - Donut style

Contain existing soft tissue around calcaneus and relieve direct pressure in heel spur, bruised and sensitive area. Available in sizes Small, Medium and Large, with 2 per package.



Extended Heel Cushions

Absorb heel strike shock and provide complete mid-foot to posterior foot cushion. Are useful in treating plantar fasciitis, heel strike impact reduction and Achilles tendonitis. Available in sizes Small, Medium and Large, with 2 per package.



Heel Wedges

Provide 15 degree declination angle. Useful for leg length discrepancy, Achilles tendonitis and heel strike impact reduction. Available in sizes Small, Medium and Large, with 2 per package.



Longitudinal Arch Supports

Support and maintain normal longitudinal arch. May be used alone in shoe or with SoftStride™ Thin Insoles. Available in sizes Small, Medium and Large, with 2 per package.



Multi-Pads w/Posting Wedges

Provide an anatomical support pad that can be used as medial/lateral posting orthotic or neutral heel cup. Multi-Pad can be modified through use of accessory Posting Wedges to provide valgus buildup. Available in sizes Small, Medium and Large, with 2 per package.



WHAT IS THE 1 MILLION STEP GUARANTEE?

Your Soft Stride insert is made from a 100% transparent medical-grade polymer with absolutely zero fillers. That means it will provide the same shock-absorbing support after your one millionth step as it does after step 1. If you notice any degradation of quality or the product has worn down, simply return it for a free replacement.

SoftStride™ Size Charts

Bone Spur Pads - Cobra style

Order No.	Size	Men	Women	EUR
71301	Small	4 - 6	5 - 7	36 - 38.5
71302	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71303	Large	11+	12+	44+

Bone Spur Pads - Donut style

Order No.	Size	Men	Women	EUR
71501	Small	4 - 6	5 - 7	36 - 38.5
71502	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71503	Large	11+	12+	44+

Extended Heel Cushion

Order No.	Size	Men	Women	EUR
71311	Small	4 - 6	5 - 7	36 - 38.5
71312	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71313	Large	11+	12+	44+

Heel Wedge

Order No.	Size	Men	Women	EUR
71341	Small	4 - 6	5 - 7	36 - 38.5
71342	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71343	Large	11+	12+	44+

Longitudinal Arch Support

Order No.	Size	Men	Women	EUR
71356	Small	4 - 6	5 - 7	36 - 38.5
71357	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71358	Large	11+	12+	44+

Multi Pad

Order No.	Size	Men	Women	EUR
71371	Small	4 - 6	5 - 7	36 - 38.5
71372	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71373	Large	11+	12+	44+

Metatarsal Pads

Order No.	Size	Men	Women	EUR
71366	Small	4 - 6	5 - 7	36 - 38.5
71367	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71368	Large	11+	12+	44+

Tenderfoot Pad

Order No.	Size	Men	Women	EUR
71411	Small	4 - 6	5 - 7	36 - 38.5
71412	Medium	6½ - 10½	7½ - 11½	39 - 43.5
71413	Large	11+	12+	44+

Thin Insoles

Order No.	Size	Men	Women	EUR
71421	A	3 - 5	4 - 6	34.5 - 37
71422	B	6 - 8	7 - 9	38.5 - 40.5
71423	C	9 - 11	10 - 12	41.5 - 44
71424	D	12 - 14		45.5 - 48
71425	E	15 - 17		48+

Sheet Products

Order No.	Dimension	Thickness
71601	14 x 9.5" / 35.5 x 24.1cm	1/10"
71604	14 x 9.5" / 35.5 x 24.1cm	1/8"

FOREFOOT PRODUCTS



Metatarsal Pads

Support transverse arch and cushion metatarsal heads. Use for generalized forefoot pain and metatarsalgia. Available in sizes Small, Medium and Large, 2 with 2 per package.



Tender Foot Pads

Separate metatarsal shafts and cushion the metatarsal heads. Unique design alleviates pain and discomfort in forefoot. Aid in treatment of neuromas and metatarsalgia. Available in sizes Small, Medium and Large, with 2 per package.

INSOLE/SHEET PRODUCTS



Thin Insoles

Reduce shock over entire plantar surface with lightweight and thinner forefoot design. Contain extra polymer in heel and arch areas to reduce shock. Available in sizes A, B, C, D or E, with 2 per package.



Sheet Products

Prefabricated of viscoelastic polymer, SoftStride™ Sheets provide an alternative to preformed foot orthotics. Useful for customization or when making an orthotic is desired. Diamond design allows for polymer flow when pressure is applied creating a softer feel. Use for prosthetic applications to provide contact cushioning against stumps, for creating customized orthotics and for shock absorbing mats where needed.

info@brownmed.com www.brownmed.com



SoftStride products are designed with the environment in mind. They are not a disposable product, therefore your purchase will reduce landfill waste over time. This product is for single patient use. Brownmed warrants this product to be free of manufacturing defects.

© 2019 Brownmed - SoftStride is a trademark of Brownmed (877) 853-5518.



MADE IN THE USA

Executive:

101 Federal Street
29th Floor
Boston, MA 02110
P 857-317-3354
F 857-317-3355

Sales:

4435 Main Street, Suite 820
Kansas City, MO 64111
TF 877-853-5518
P 816-581-7001
F 816-581-7005

Manufacturing:

1300 Lundberg Drive West
Spirit Lake, IA 51360
P 712-336-4395
F 712-336-2874