



FOOT/ANKLE WRAP

UNIVERSAL + ADJUSTABLE



CONFORMS WELL TO FOOT/ANKLE



PAIN RELIEF FOR BRUISES, SPRAINS OR INJURIES REQUIRING COLD THERAPY



DELIVERS BENEFITS OF CRYOTHERAPY



PROVIDES COMPRESSION FOR SOFT TISSUE TRAUMA



SOFT FLEECE PROTECTS SKIN



FOOT/ANKLE WITH ICE PACK INSERT



Polar Ice Foot/Ankle Wrap delivers the benefits of cryotherapy with compression for soft tissue trauma. By reducing pain and controlling edema or swelling, cold therapy allows faster return of normal function. Encapsulated ice packs filled with purified water are encased in soft fleece material that eliminates migration of ice wrap during use. Easily adjusted with hook & loop closure straps, fleece wrap conforms to body and insulates skin to prevent frost bite. Use post-surgically, during rehabilitation and following physical activity.

REVERSIBLE TO FIT RIGHT OR LEFT FOOT/ANKLE

30405 UNIVERSAL SIZE



www.brownmed.com



Executive

101 Federal Street, 29th Floor
Boston, MA 02110
P 857-317-3354
F 857-317-3355

Sales

4435 Main Street, Suite 820
Kansas City, MO 64111
TF 877-853-5518
P 816-581-7001
F 816-581-7005

Manufacturing

1300 Lundberg Drive West
Spirit Lake, IA 51360
P 712-336-4395
F 712-336-2874

2019-0100 FAW