

# Herbst Cradle<sup>™</sup>



Patent No. 5,876,364

## Unique Design

- lightweight, cool and soothing to irritated skin
- antibacterial filling
- cradle alleviates pressure to the heels and malleoli
- anti-rotation bar keeps the foot at a 90 degree angle
- strategically placed bead pockets stimulate blood circulation
- hook and loop fasteners allow for secure and expedient application

## Used For

- prevention and healing of decubital ulcers on the heel and malleoli
- prevention of bed linen abrasions
- medial or lateral rotations of the lower extremities and unwanted hip positioning
- plantar flexion contractures (Achilles' contractures)

Item #	Size	Women	Men
04020	Small	≤6	≤5
04021	Medium	6-10.5	5-9.5
04022	Large	>10.5	>9.5
50043	X-Large		

### Executive:

101 Federal Street  
29th Floor  
Boston, MA 02110  
P 857-317-3354  
F 857-317-3355

### Sales:

4435 Main Street, Suite 820  
Kansas City, MO 64111  
TF 877-853-5518  
P 816-581-7001  
F 816-581-7005

### Manufacturing:

1300 Lundberg Drive West  
Spirit Lake, IA 51360  
P 712-336-4395  
F 712-336-2874